

Creatinase (C2-AE)

from recombinant *E. coli*

Creatine amidinohydrolase, EC 3.5.3.3



SPECIFICATION


Appearance	white lyophilizate
Activity	≥9 U/mg lyophilizate
Contaminant	catalase ≤0.5%
Stabilizer	sucrose
Storage	at -20°C

PROPERTIES

Molecular weight	ca. 80 kDa (gel filtration)
Structure	2 subunits of 46 kDa (SDS-PAGE)
Michaelis constant	1.3×10^{-2} M (creatine)
pH Optimum	7.0–9.0 (Fig. 1)
pH Stability	5.0–11.0 (Fig. 2)
Optimum temperature	40°C (Fig. 3)
Thermal stability	below 45°C (Fig. 4)
Stability (liquid form)	stable at 37°C for at least two weeks (Fig. 5)
Stability (powder form)	stable at 30°C for at least one month (Fig. 6)
Inhibitor	Hg ²⁺

ASSAY PROCEDURE

Principle



The appearance of urea is measured spectrophotometrically at 435 nm.

Definition of unit

One unit (U) is defined as the amount of enzyme which produces 1 μmol of urea per min at 37°C and pH 7.7 under the conditions described below.

Reagents

- Potassium phosphate buffer, 0.3 M; pH 7.7: mix 0.3 M KH_2PO_4 solution and 0.3 M K_2HPO_4 solution to make a pH 7.7 solution.
- Creatine solution, 0.1 M: dissolve 1.49 g of creatine monohydrate/100 ml of distilled water.
- p*-Dimethylaminobenzaldehyde (DMAB) solution, 0.87%: dissolve 2.0 g of DMAB in 100 ml of ethanol (99.5%) and add 15 ml of conc. HCl and 115 ml of distilled water.
- Enzyme dilution buffer: mix 10 mM KH_2PO_4 solution and 10 mM K_2HPO_4 solution to make a pH 8.0 solution and add 2-mercaptoethanol (0.16 ml/l of the buffer).

Sample: dissolve the lyophilized enzyme to a volume activity of 1.4–2.8 U/ml with ice-cold enzyme dilution buffer (Reagent D) immediately before measurement.

Procedure

- Pipette the following reagents into a test tube.

0.1 ml	Potassium phosphate buffer	(Reagent A)
0.8 ml	Creatine solution	(Reagent B)
- Equilibrate at 37°C for about 5 min.
- Add 0.1 ml of sample and incubate for 10 min at 37°C.
- Add 2 ml of DMAB solution (Reagent C) and allow to stand for about 30 min at 25°C.
- Read the absorbance at 435 nm in a cuvette (light path: 1 cm) (A_s).

The blank solution is prepared by reversing the sequence of addition of sample and DMAB solution (Reagent C) (A_0).

Calculation

Activity can be calculated by using the following formula:

$$\text{Volume activity (U/ml)} = \frac{(A_s - A_0) \times 1.0(\text{ml}) \times df}{0.0543 \times 0.1(\text{ml}) \times 10(\text{min})} = \Delta A \times 18.4 \times df$$

$$\text{Weight activity (U/mg)} = (\text{U/ml}) \times 1/C$$





0.0543 : Millimolar extinction coefficient of urea under the assay conditions ($\text{cm}^2/\mu\text{mol}$)

df : Dilution factor

C : Content of creatinase preparation in sample (mg/ml)

APPLICATIONS

The enzyme is useful for the determination of creatinine and creatine in clinical analysis.

**REFERENCES**

Suzuki, M, *Medical Technology*, **7**, 945–950 (1979).

Suzuki, M. and Yoshida, M., *Clin. Chim. Acta*, **140**, 289–294 (1984).

Suzuki, M. and Yoshida, M., *Clin. Chim. Acta*, **143**, 147–155 (1984).

EXPERIMENTAL DATA

Fig. 1 pH Optimum




Fig. 2 pH Stability




Fig. 3 Optimum temperature




Fig. 4 Thermal stability




Fig. 5 Stability (liquid form) at 37°C




Fig. 6 Stability (powder form) at 30°C

