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Product code:60337

**CAUTION !**

1. Do not attempt to drink the kit reagent or touch it with bare hands or allow it to splash into eyes.
2. Please make sure to read the precautions and instructions in this Instruction Manual before attempting to use the kit and exercise extreme caution when using it.

LuciPac LS2.8-400 is a kit for testing cleanliness levels mainly for canalicular objects using bioluminescence techniques using firefly luciferase developed with Kikkoman's unique biotechnology.

This kit contains long stem swabs suitable for swabbing inside of the channels of medical endoscopes and of filling nozzles in beverage industry, which have not been accessible by LuciPac Pen and/or other ATP swab test kits.

Size of the long stem swab: Cotton bud 2.8mmΦ and stem length 400 mm

[Applications]

This kit can be used to test cleanliness levels using a swab test. Do not use this device for any purpose other than what it was designed for as a tester of cleanliness levels. Please be aware that this kit cannot be used to test or measure for counts of viable bacteria or more specific types of pathogenic bacteria.

This kit is used to indicate cleanliness levels in food production centers and medical treatment facilities. This kit shall not provide any guarantee that a given test environment is free of bacterial contamination. This kit cannot be used to test food or food products directly. This kit is not a medical device.

[Measurement Principles]

This kit uses an enzyme cycling method based on a combination of luminescent reactions from firefly luciferase and pyruvate, orthophosphate dikinase (PPDK). This method produces a given amount of luminescent that is proportional to the amounts of adenosine triphosphate (ATP) and adenosine monophosphate (AMP) present.


Fig. 1. Principle of luminescence method using the enzyme cycling

Firefly luciferase emits light in the presence of ATP and luciferin. The AMP produced from this reaction is converted back into ATP using PPDK to enable a high but stable amount of luminescence to be obtained.

ATP is a source of energy necessary for various forms of life that is present in organic residues, such as microorganisms, food residue, and biological substances that originate from other living organisms. This method allows you to measure and detect organic

residues quickly and at high sensitivity by detecting ATP using luciferase, which is why it is widely used in determining cleanliness levels in food production centers and medical treatment facilities.

In addition, this kit can be used to measure amounts of not only ATP but also of AMP, which is produced from the breakdown of ATP, enabling it to be used to test for an even wider range of organic residues.

[Contents]

The LuciPac LS2.8-400 kit comes with 40 individually packed long stem swabs and two aluminum bags each containing 20 LuciPac Pen-AQUA devices (for a total of 40 devices).



Fig. 2. Each part name of the long stem swab and the LuciPac Pen-AQUA.

Table 1. Main components of each reagent

Reagent name	Main component
Luminescent reagent	Luciferin Luciferase Magnesium acetate Phosphoenolpyruvic acid Pyrophosphoric acid Pyruvate,orthophosphate dikinase
Releasing reagent	Surfactant (Benzalkonium chloride)

[Precautions for Use]

Please make sure to follow the items outlined below in order to obtain optimal performance from this kit.

- ① Please take out two aluminum bags of “LuciPac Pen-AQUA” immediately upon receive this box and keep in a refrigerator (2-8 °C). The long stem swabs may be kept at room temperature.
- ② The long stem swab was specially manufactured for testing cleanliness. Other commercial cotton swab may not produce correct results.
- ③ Choose LuciPac LS2.8-400 or LS3.2-400, depending on the inner diameter of the object to be tested.
- ④ Do not push the long stem swab into the narrow inspection point forcibly. Stop pushing it forcibly when there is strong resistance at the time of the insertion. If not done, the cotton bud may fall off when you pull the swab up.
- ⑤ Do not use products that are past the expiry date. Expired products may not yield accurate results (Expiration date is printed on the labels of the kit box and aluminum bags.)
- ⑥ Be sure to use only designated products when taking measurements for amounts of luminescence. This kit cannot be used with nondesignated products.
- ⑦ The LuciPac Pen-AQUA should be allowed to reach room temperature (see Table 2) prior to testing, if they are from refrigerated stock. Measurement values may read lower than actual if the LuciPac Pen-AQUA devices are used while still cold. Use the LuciPac Pen-AQUA as soon as possible once they have returned to room temperature. Do not let the kit sit out in temperatures exceeding 35°C (95°F). High temperatures may cause product performance to drop.

Table 2. Appropriate measurement temperature

Model	Temperature compensation setting	Temperature range
Lumitester PD-20	—	20 - 35°C (68 - 95°F)
Lumitester PD-30	OFF	
		ON*1)

*1) Temperature compensation is performed in accordance with the measured value of a thermometer incorporated in the main body of Lumitester PD-30. Consequently, the compensation cannot be performed precisely if the temperature of the main body of PD-30 is different from LuciPac Pen-AQUA. Please be sure that temperature of the main body of PD-30 and LuciPac Pen-AQUA become same before use.

- ⑧ Wherever possible, be sure to use up all LuciPac Pen-AQUA devices from a single bag that has been opened at one time. If you absolutely must store leftover LuciPac Pen-AQUA devices once finished with a test session, firmly close the aluminum bag and store it in a refrigerated environment (2°C to 8°C (35.6°F to 46.4°F)). High temperatures may cause product performance to drop.
- ⑨ Do not subject the kit or any part of it to direct sunlight for long periods of time. Strong light may cause product performance to drop.
- ⑩ Do not touch any of the parts inside the LuciPac Pen-AQUA and tip of the long stem swab, particularly not cotton bud of

the long stem swab and any part of the sampling stick (red part), with fingers or other object before use. Touching the parts may affect cleanliness levels, making them hard to determine.

- ⑪ Do not drop the kit or any of its parts or allow any parts to be struck or jolted. The inner aluminum sheets and other parts in the kit may become damaged, causing product performance to drop.
- ⑫ Do not use the kit if any parts become damaged such as the inner aluminum sheets. Such damage may affect product performance, causing cleanliness levels to fail to be measured correctly. You can tell if the aluminum sheet is damaged or not by checking to see if the releasing reagent is leaking.
- ⑬ Secure the test tube of the LuciPac Pen-AQUA to ensure that it does not become shifted or displaced. It may be difficult to accurately determine cleanliness levels if the test reagent begins leaking because the tube has been shifted. It may become difficult to remove the LuciPac Pen-AQUA from the measurement device once the measurements have been taken. Moreover, it may cause malfunction of measurement device.

[Designated Product for Measurement]

Lumitester PD-20 and Lumitester PD-30

(Manufacturer: Kikkoman Biochemifa Company)

Be sure to use only designated products when taking measurements.

[Measurement Methods]

1. Measurement procedures

Complete the procedures listed below within appropriate temperature range as shown in Table 2 depending on an instrument used and temperature compensation setting. Make sure to always run measurement tests at the same temperature to maintain repeatability from comparison, if the temperature compensation is not used. Remove the LuciPac Pen-AQUA from refrigerator, and wait until they have reached room temperature (see Table 2). Use them as soon as possible once they have returned to room temperature.

- ① Wear lab gloves to prevent contamination and take out the long stem swab carefully so as to not to touch the parts as this may affect the cleanliness level from the end of the stem of the individual package.
- ② If running the measurement test on a dry surface, first moisten the cotton bud or the test subject with tap water*2). Please avoid using saline solution instead of water.
- ③ Handle the long stem swab with care so as to not to touch the parts as this may affect the cleanliness level. Insert the long stem swab in the canalicular object to be tested, and wipe it off as long as the length of the long stem swab reaches, and pull the swab up slowly.
- ④ Hold the long stem swab at 12 cm or 13 cm distant point from the cotton bud, and remove the sampling stick from the main body (casing) of the LuciPac Pen-AQUA by your free fingers. Handle the sampling stick with care so as to not to touch the parts as this may affect the cleanliness level.
- ⑤ Insert the long stem swab in the main body (casing) of the LuciPac Pen-AQUA to push and break the center of the upper seal of the container for releasing reagent. Be careful not to break bottom seal of the container for releasing reagent. Wash the cotton bud in the releasing reagent. Remove the long stem swab after thoroughly washing the cotton bud in the releasing reagent.
- ⑥ Return the sampling stick to the main body (casing) and push it all the way into the main body (casing).
- ⑦ Hold firmly onto the LuciPac Pen-AQUA casing and shake it. Allow leftover luminescent reagent to thoroughly dissolve.
- ⑧ Insert the LuciPac Pen-AQUA into the Lumitester and close the Lumitester measurement chamber cover.

© 2014 Kikkoman Corp. (20150702)

- ⑨ Press ENTER key. Measurement result is displayed after countdown.

*2) Regular tap water may contain slight amount of ATP and AMP. We recommend that you wash or otherwise clean off the faucet in such cases and allow the water to run for a while before collecting any for use.

2. Handling of data

Normal / defect criteria should be decided by user based on observation of data for certain period of time under normal operation in consideration with fluctuation of the data. In addition, if possible, artificial inadequate cleaning conditions or artificially contaminated conditions should be evaluated.

[Disposal Methods]

This kit doesn't contain any GHS hazardous chemical. This kit can be disposed of as regular garbage, but when disposing of it, it would be better to separate the parts and dispose of each one properly in accordance with the local regulations outlined by the local governments for proper disposal of waste materials. Moreover, when the kit was used to test medical devices, dispose it as clinical waste if needed.

The main materials and parts used in this kit are listed below. No PVC materials are used in the making of the plastics in this kit.

Table 3. Main raw materials of the structural parts of this product

	Structural parts	Raw materials
LuciPac Pen-AQUA	Stick holder (orange part)	Polypropylene
	Main body (casing)	Polypropylene
	Sampling stick (red part)	ABS resin
	Container for releasing reagent	Polypropylene, Aluminum
	Measurement tube	Polypropylene, Aluminum
Long stem swab	Aluminum bag (with dehumidifying function)	Aluminum, Polyethylene, Polyethylene terephthalate
	Stem	Polyethylene
	Cotton bud	Cotton
	individual package	Paper, Polyester, Polypropylene
External packaging	Outer package	Paper
	Shrink-wrap	Polypropylene

[Precautions for Handling]

- ① Be careful not to let the reagents or other substances in the kit get into your mouth or eyes, or onto bare hands before or after use. Rinse your mouth out thoroughly with water if any of the substances get into your mouth, rinse off your skin with copious amounts of water if any get onto your skin, and rinse out affected eyes thoroughly with copious amounts of water should it get into your eyes, and immediately contact a physician for advice and follow the instructions given.
- ② Exercise enough caution when storing and disposing of the kit and its reagents to ensure that none of the substances become mixed in with food and other products.
- ③ Be careful not to get fingers caught when inserting the sampling stick into the main body (casing).
- ④ Please make sure to store this kit and its parts out of the reach of young children.
- ⑤ Note that the releasing reagent used in this kit contains cationic surfactant (benzalkonium chloride^{*3)}). Take precautions when disposing of this kit after use to ensure that such substances do not become mixed in with food products at food production centers and similar facilities.

*3) Benzalkonium chloride is a disinfectant and antiseptic commonly used in hand and finger sterilizer solutions.

[Storage]

- ① Kit storage: The aluminum bags of LuciPac Pen-AQUA are to be stored at low temperature (2°C to 8°C (35.6°F to 46.4°F)) for long term storage. The LuciPac Pen-AQUA can be stored below 25°C (77°F) for up to 14 days or below 30°C (86°F) for up to 5 days before opening an aluminum bag without any adverse effect on the long term stability. Do not freeze the LuciPac Pen-AQUA. The long stem swabs may be kept at room temperature.
- ② We recommend that you use all 20 LuciPac Pen-AQUA devices in a single bag at one time after opening an aluminum bag. If you have leftover LuciPac Pen-AQUA devices that you must store after opening a bag, be sure to store them at the recommended low temperature (2°C to 8°C (35.6°F to 46.4°F)) and use them within two weeks from when the bag is opened.
- ③ Expiry date: Printed on the labels of the outer box and aluminum bag.

[Warranty]

Kikkoman Biochemifa Company warrants the products in this kit to have a certain level of quality. This warranty guarantees that Kikkoman Biochemifa Company shall replace defective products should any be found. This warranty does not provide any other guarantees. Kikkoman Biochemifa Company shall not be liable for any damages, including special or consequential damages, or expenses arising directly or indirectly from the use of this product.

Symbols used in the packaging and labeling of this product	
	Symbol for "temperature limitation." The upper and lower temperature limits will be indicated on either side of the symbol. Please store this product at the indicated temperature range.
	Symbol for "Caution" or "Attention" for use.
	Symbol for "Lot Number." This symbol shall be adjacent to the manufacturer's lot number (e.g. 20150410Q) or description of its printed location.
	Symbol for "Use By." This symbol shall be adjacent to the expiration date, expressed as YYYYMMDD (e.g. 20160709), or description of its printed location.
	Symbol for "Consult Instructions Manual."
	Symbol for "Manufacturer." This symbol shall be adjacent to the name and address of the manufacturer.

"LuciPac" and "Lumitester" are registered trademarks of Kikkoman Corporation in Japan and other countries.

Manufacturer:


Kikkoman Biochemifa Company

2-1-1 Nishi-shinbashi, Minato-ku, Tokyo 105-0003, Japan
 Phone: +81-3-5521-5490 / Fax: +81-3-5521-5498
 E-mail: biochemifa@mail.kikkoman.co.jp
 URL: <http://biochemifa.kikkoman.co.jp/e/>


© 2014 Kikkoman Corp. (20150702)

[Measurement Methods]


Leave "LuciPac Pen-AQUA" at room temperature until they have reached room temperature.


Wear lab gloves to prevent contamination and take out the long stem swab carefully so as to not to touch the parts as this may affect the cleanliness level.


If running the measurement test on a dry surface, first moisten the cotton bud or the test subject with tap water*4). Please avoid using saline solution instead of water.


Handle the long stem swab with care so as to not to touch the parts as this may affect the cleanliness level. Insert the long stem swab in the canalicular object to be tested, and wipe it off as long as the length of the long stem swab reaches, and pull the swab up slowly.


Hold the long stem swab at 12 cm or 13 cm distant point from the cotton bud, and remove the sampling stick from the main body (casing) of the LuciPac Pen-AQUA by your free fingers. Handle the sampling stick with care so as to not to touch the parts as this may affect the cleanliness level.

Insert the long stem swab in the main body (casing) of the LuciPac Pen-AQUA to push and break the center of the upper seal of the container for releasing reagent. Remove the long stem swab after thoroughly washing the cotton bud in the releasing reagent.

Return the sampling stick to the main body (casing) and push it all the way into the main body.

Hold firmly onto the LuciPac Pen-AQUA casing and shake it. Allow leftover luminescent reagent to thoroughly dissolve.

Insert the LuciPac Pen-AQUA into the Lumitester and close the Lumitester measurement chamber cover.

Press ENTER key. Measurement result is displayed after countdown.

*4) If the luminescence level of tap water shows high, we recommend you to let the water flow a few minutes and/or clean up the outlet of tap.