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ATP Hygiene Monitoring test kit (Swab Test)

LuciPac II Instruction Manual

Product code:61301

CAUTION!

- 1. Do not attempt to drink the kit reagent or touch it with bare hands or allow it to splash into eyes.
- 2. Please make sure to read the precautions and instructions in this Instruction Manual before attempting to use the kit and exercise extreme caution when using it.

The Kikkoman ATP hygiene monitoring test kit LuciPac II is a kit for testing cleanliness levels using bioluminescence techniques using firefly luciferase developed with Kikkoman's unique biotechnology.

[Applications]


This kit can be used to test cleanliness levels using a swab test.

Do not use this device for any purpose other than what it was designed for as a tester of cleanliness levels. Please be aware that this kit cannot be used to test or measure for counts of viable bacteria or more specific types of pathogenic bacteria.

This kit is used to indicate cleanliness levels in food production centers and medical treatment facilities. This kit shall not provide any guarantee that a given test environment is free of bacterial contamination. This kit cannot be used to test food or food products directly.

[Measurement Principles]

This kit uses a luminescent reactions from firefly luciferase.


Firefly luciferase emits light in the presence of ATP and luciferin.

ATP is a source of energy necessary for various forms of life that is present in organic residues, such as microorganisms, food residue, and biological substances that originate from other living organisms. This method allows you to measure and detect organic residues at high speed and high sensitivity by detecting ATP using luciferase, which is why it is widely used in determining cleanliness levels in food production centers and medical treatment facilities.

This kit is a simple integrated testing instrument that contains both the test reagent and the swab device required for testing cleanliness levels.

[Contents]

The LuciPac II kit comes with five aluminum bags each containing 20 swab devices (for a total of 100 swab devices).

Reagent name	Main component
Luminescent reagent	Luciferin Luciferase Magnesium acetate
Releasing reagent	Surfactant (Benzalkonium chloride)

[Precautions for Use]

Please make sure to follow the items outlined below in order to obtain optimal performance from this kit.

- ① Do not use products that are past the expiry date. Expired products may not yield accurate results (the expiry date is printed on the bottom of the aluminum bag holding the swab devices).
- ② Be sure to use only designated products when taking measurements for amounts of luminescence. This kit cannot be used with nondesignated products.
- 3 The swab devices should be allowed about 20min to reach room temperature(20°C to 35°C (68°F to 95°F)) prior to testing, if they are from refrigerated stock. Measurement values may read lower than actual if the swab devices are used while still cold. Use the swab devices as soon as possible once they have returned to room temperature. Do not let the kit sit out in temperatures exceeding 35°C (95°F). High temperatures may cause product performance to drop.
- Wherever possible, be sure to use up all swab devices from a single bag that has been opened at one time. If you absolutely must store leftover swab devices once finished with a test session, firmly close the aluminum bag and store it in a refrigerated environment (2°C to 8°C (35.6°F to 46.4°F)). High temperatures may cause product performance to drop.
- Do not subject the kit or any part of it to direct sunlight for long periods of time. Strong light may cause product performance to drop.
- 6 Do not touch any of the parts inside the kit, particularly not any part of the cotton swab itself, with a finger or other object before use. Touching the parts may affect cleanliness levels, making them hard to determine.
- ⑦ Do not drop the kit or any of its parts or allow any parts to be struck or jolted. The inner aluminum sheets and other parts in the kit may become damaged, causing product performance to drop.
- 8 Do not use the kit if any parts become damaged such as the inner aluminum sheet. Such damage may affect product performance, causing cleanliness levels to fail to be measured correctly. You can tell if the aluminum sheet is damaged or not by checking to see if the releasing reagent is leaking.
- Set the benchmark level at which the required cleanliness levels can
 be obtained. If the benchmark level is not set correctly, cleanliness
 levels may not be evaluated correctly.
- ⑤ Secure the test tube of the LuciPac II ensure that it does not become shifted or displaced before measurement. It may be difficult to accurately determine cleanliness levels if the test reagent begins leaking because the tube has been shifted.

[Designated Product for Measurement]

Lumitester C-110 and Lumitester C-100

(Manufacturer: Kikkoman Biochemifa Company)

Be sure to use only designated products when taking measurements.

[Measurement Methods]

1. Measurement procedures

Complete the procedures listed below within a temperature range of 20°C to 35°C (68°F to 95°F). Make sure to always run measurement tests at the same temperature to maintain repeatability from comparison. Remove the LuciPac II from refrigerator, and wait until they have reached room temperature (about $20~\text{min}, 20^{\circ}\text{C}$ to 35°C (68°F to 95°F)). Use the swab devices as soon as possible once they have returned to room temperature.

- ① Remove the swab stick from the main body (casing).
- ② If running the measurement test on a dry subject, first moisten the cotton swab or the test subject with tap water *1).
- ③ Use it to swab the test subject with constant pressure as much as possible. Under suitable pressure, the swab shaft will slightly bend as illustrated in the panel 3 of page 6.
- 4 Return the swab stick to the the main body (casing) and push it all the way into the the main body (casing).
- (5) Hold firmly onto the LuciPac II casing and shake it.
- 6 Allow the leftover luminescent reagent to thoroughly dissolve.
- The put your thumb on one of the two welded zones and pull the tip of the test tube to break the welded bound. Repeat this procedure for another welded zone of the opposite side.
- 8 Remove the test tube from the main body. Be careful not to spill the releasing reagent.
- ⑤ Insert the test tube into the measuring chamber of Lumitester. Press the "ENTER" key. Results are obtained in 10 seconds.
- *1) Regular tap water may contain slight amount of ATP. We recommend that you wash or otherwise clean off the faucet in such cases and allow the water to run for a while before collecting any for use.

Refer to "Direction for LuciPac II" on page 6 for further details.

2. Handling of data

Benchmark levels 1, 2 are set that are used to determine whether cleanliness levels are at acceptable levels or not (see Fig. 1). Cleanliness levels are designated as Pass (A) if the measured amount of luminescence falls below level 1 and as Fail (C) if the amount is higher than level 2. With Fail (C), cleaning will need to be redone and cleaning procedures may need to be revised. Cleanliness levels where the measured amount of luminescence is above level 1 but below level 2 are designated as Caution (B), and it is recommended that cleaning be redone and cleaning procedures revised in such cases.

Examples of benchmark setting values for food production centers are listed in Table 1. The benchmark setting values for cleanliness level control shown in Table 1 are provided merely for reference and may not necessarily apply to all situations and cases. The benchmark levels are important guidelines for determining cleanliness levels and should be set as appropriate according to required cleanliness levels.


Fig. 1: Setting Benchmark Level for Cleanliness Level Control

Table 1: Example of Setting Benchmark Level for Cleanliness Level Control in Food Production Centers

Relative Light Unit (RLU)	
Level 1	Level 2
100	200
500	1,000
1,500	3,000
	Level 1 100 500

[Disposal Methods]

This kit contains no hazardous materials. This kit can be disposed of as regular garbage, but when disposing of it, it would be better to separate the parts and dispose of each one properly in accordance with the local regulations outlined by prefecture and city, town, or village for proper disposal of waste materials.

The main materials and parts used in this kit are listed below. No PVC materials were used in the making of the plastics in this kit.

Structural parts	Raw materials		
Swab stick (red part)	polypropylene		
Main body (casing)	polypropylene		
Swab shaft	polypropylene		
Container for releasing reagent	polypropylene		
Measurement tube	polypropylene		
Desiccant	Silica gel, polyethylene, polypropylene		
Aluminum bag	Aluminum, polyethylene, polyethylene terephthalate		
Outer bag	Polyethylene		

[Precautions for Handling]

Please observe the following items to ensure safe use of this product.

- ① Be careful not to let the reagents or other substances in the kit get into your mouth or eyes, or onto bare hands before or after use. Rinse your mouth out thoroughly with water if any of the substances get into your mouth, rinse off your skin with copious amounts of water if any get onto your skin, and rinse out affected eyes thoroughly with copious amounts of water should it get into your eyes, and immediately contact a physician for advice and follow the instructions given.
- ② Exercise enough caution when storing and disposing of the kit and its reagents to ensure that none of the substances become mixed in with food and other products.
- 3 Be careful not to get fingers caught when inserting the swab stick into the main body (casing).
- Please make sure to store this kit and its parts out of the reach of young children.
- S Note that the releasing reagent used in this kit contains cationic surfactants (benzalkonium chloride*2)). Take precautions when disposing of this kit after use to ensure that such substances do not become mixed in with food products at food production centers and similar facilities.
- *2) Benzalkonium chloride is a disinfectant and antiseptic commonly used in hand and finger sterilizer solutions.

[Storage]

- 1) Kit storage: Kits are to be stored at a low temperature (2°C to 8°C (35.6°F to 46.4°F)) for long term storage. The kit can be stored below 25°C (77°F) for up to 14 days or below 30°C (86°F) for up to 5days before opening an aluminum bag without any adverse effect on the long term stability. Do not freeze the kit.
- 2) We recommend that you use all 20 swab devices in a single bag at one time after opening an aluminum bag. If you have leftover swab devices that you must store after opening a bag, be sure to store them at the recommended low temperature (2°C to 8°C (35.6°F to 46.4°F)) and use them within two weeks from when the bag is opened.
- 3) Expiry date: Printed on the aluminum bag.

[Warranty]

Kikkoman Biochemifa Company warrants the products in this kit to have a certain level of quality. This warranty guarantees that Kikkoman Biochemifa Company shall replace defective products should any be found. This warranty does not provide any other guarantees. Kikkoman Biochemifa Company shall not be liable for any damages, including special or consequential damages, or expenses arising directly or indirectly from the use of this product.

Manufacturer:


Kikkoman Biochemifa Company

2-1-1 Nishi-shinbashi, Minato-ku, Tokyo 105-0003, Japan Phone: +81-3-5521-5490 / Fax: +81-3-5521-5498

E-mail: biochemifa@mail.kikkoman.co.jp URL: http://biochemifa.kikkoman.co.jp/e/


[Direction for LuciPac II]

<Leave "LuciPac II" at room temperature (20-35°C) for 20 minutes before use >


Pull the cotton swab stick
(A) out of the main body


When the test object is dry, moisten the swab (A) with tap water* or wet the object with water*

Swab the test object with the swab stick (A).


Put the swab stick (A) back to the main body (B) and push it through all the way by putting the tip of the test tube (C) on a palm of hand or table to avoid having the tube (the tip part of the body) come off the main body.


Shake the whole body a few times so that any remaining liquid in the capsule falls into the test tube (C). Otherwise, the remaining liquid may spill out when you remove the test tube in step 8.


Gently shake the whole body so that the luminescent reagent is entirely dissolved. Otherwise, the resulting value may become lower than the accurate value.

Put your thumb on one of the two welded zones and pull the tip of the test tube to break the welded bound. Repeat this procedure for another welded zone of the opposite side.

Remove the test tube from the main body. Be careful not to spill the releasing reagent.

Insert the test tube into the measuring chamber of Lumitester. Press the "ENTER" key. Results are obtained in 10 seconds.

If the luminescence level of tap water shows high, we recommend you to let the water flow a few minutes and/or clean up the outlet of the tap.